

NOBAMED Paul Danz AG

SCHUTZUNTERLAGE

REF 670079

Product Description, Intended use, Application

The protection sheet, **size 75 x 90 cm**, is provided with a green plastic coating which is impermeable to liquids. The absorbent protection sheet rapidly absorbs liquids, while the foil coating on its underside prevents that the sheet is soaked. The protection sheet is used in medical sectors and in nursing to protect stretchers, couches, beds, or tables from soiling, thus reducing the risk of contamination and counteracting the spread of germs. The protection sheet is for single use.

Composition

Cellulose, polyethylene

Contraindications

No contra-indications known.

The product should not be used if there is a known hypersensitivity to one of the belowmentioned ingredients.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of

cellulose

Symbols used in labelling

Explanations at www.nobamed.com







Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN 1041.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: March 12, 2021 [REV 5] replaced: April 01, 2020 [Rev 4]