



NOBAMED Paul Danz AG

NOBAPROTECT® FACE SHIELD

REF 672090 (from LOT 137397)

Product Description and Purpose

Face shield in accordance with EN 166:2001, size: approx. 30 cm x 20 cm, thickness: 0,5 mm. The field of application of the product is the protection against droplets and splashes of liquids in healthcare. It complies with the optical requirements for face shields, is without sharp edges and impairments of the visual field, is resistant to fogging, non-flammable and does not cause skin irritations. The face shield is suitable for spectacle wearers and can be worn over a visual aid. Similarly, a surgical or FFP mask can be worn under the face shield. When used with other protective equipment (hearing protection, respiratory protection, etc.), the face shield must be put on in such a way that the protective effects of the other protective equipment are not impaired in their effectiveness or interfere with each other.

Composition

PET, polyurethane

Notes

Read the instructions carefully before using the product.

Normative and Regulative Requirements

NOBAPROTECT® FACE SHIELD is PPE CAT II according to the PPE Directive (EU) 2016/425 and complies with the requirements of EN 166:2001.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary: plastic bag
Secondary: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



1N-166 3S CE

PSA CAT II (EU) 2016/425

Storage and transport

If the product is wet, allow it to dry completely before storing it. To avoid scratches, make sure that the product is not lying with the viewing window facing downwards. Store it in a well-ventilated place, away from extreme temperatures and direct sunlight. Never place heavy objects on it. Avoid excessive folding and store preferably hanging.

Disposal

In accordance with applicable legislation and in compliance with infection prevention requirements