



## WASCHHANDSCHUH

REF 672200

### Product Description, Intended use, Application

The wash glove for use on both sides, is made of two soft layers of non-woven fabric. The inside of the gloves is not plasticized. As the layers are sealed together at the edge, the glove is extremely durable. The offset opening makes it easy to put on and take off. The WASCHHANDSCHUH is used to clean the patient's body to prevent the spread of germs.

### Composition

Polyester

### Contra-Indications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.