

# NOBAMED Paul Danz AG

## **WASCHHANDSCHUH-basic**

#### **REF 672201**

## **Product Description, Intended use, Application**

The WASCHHANDSCHUH-basic consists of two sealed soft and fine layers of Molton-like non-woven fabric. The inside of the wash glove is not plasticized. The offset opening makes it easy to put on and take off. The WASCHHANDSCHUH-basic is used to clean the patient's body to prevent the spread of

### Composition

Polyester

#### **Contra-Indications**

The product should not be used in the case of a known allergy against the material.

#### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **Normative and Regulative** Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

#### **Packaging**

Primary packaging: PE-bag

Secondary packaging: carton made of

cellulose

#### Symbols used in labelling Explanations at www.nobamed.com





Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

#### **Storage and Transport**

To be stored in a dry and dust-free environment

#### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.