

REF 694520 (as of LOT 137715)

Product Description, Intended use, Application

The sterile oral hygiene pouch consists of 10 x loop gauze balls, \varnothing 3,0 cm

1 x universal clamp, shear-shaped, blue

1 x plastic tray with closing lid

1 x plastic cup

The gauze balls can be soaked with a suitable solution for oral hygiene. By means of the universal clamp the patient's oral cavity can be cleaned with the balls.

Composition

Loop gauze balls: cotton Universal clamp: polystyrene Plastic tray and cup: PET

Contrainidcations

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN ISO 11135.

The product does not contain dangerous toxic substances according to REACH.

Packaging:

Primary packaging: paper-film-

packaging

Secondary packaging: PE-bag

Tertiary packaging: carton made of

cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: February 28, 2024 [Rev 6], replaces: August 03, 2022 [Rev 5]