



Dressing change set with forceps

REF 696001

Product Description, Intended use, Application

The sterile dressing change set consists of 5 gauze balls \varnothing 3 cm and 5 gauze swabs 10 x 10 cm. They are placed in a plastic tray, together with blue plastic forceps. The set is wrapped in a water resistant sheet.

Composition

Cotton, silicone, polyethylene, polystyrene, cellulose

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary: paper-film-packaging

Secondary: PE bag

Tertiary: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree, protect from sunlight.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.