



## NOBALINE<sup>®</sup>-Bettuch

REF 772580

### Product Description, Intended use, Application

NOBALINE<sup>®</sup>-Bettuch is an **80 cm x 250 cm** aluminium-coated non-woven bed sheet that is non-adhesive to wounds. The sterile bed sheet is mainly used as an underlay for the treatment of large wounds and abrasions in intensive care. The NOBALINE<sup>®</sup>-Bettuch can be used as an overlay for normal beds as well as for special beds. NOBALINE<sup>®</sup>-Bettuch is intended for single use and is individually sealed.

### Composition

Polyethylene, aluminum, polyester, rayon

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: paper-film packaging

Secondary packaging: PE-bag

Tertiary packaging: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.