

# NOBAMED Paul Danz AG

# **VERBANDTUCH-steril**

### **REF 776406**

# **Product Description; Intended use, Application**

The wound sheet, **size BR** (small) measuring **40 cm x 60 cm**, is used for aseptic wound coverage and for pressure-free coverage of large-scale injuries. VERBANDTUCH-steril is for single use, especially in first-aid, and is individually sealed in a paper-film packaging.

# Composition

Polypropylene

#### **Contra-Indications**

The product should not be used in the case of a known allergy against the material.

#### **Notes**

DIN 13152-BR 400 x 600 mm

#### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The wound sheet complies with the normative requirements of DIN 13152.

Sterilization in compliance with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: paper-film

packaging

Secondary packaging: PE bag

Tertiary packaging: carton

made of cellulose

### Symbols used in labelling

Explanation at www.nobamed.com









Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### **Storage and Transport**

Dry and dustfree, protect from sunlight.

### Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

# Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: September 13, 2022 [REV 3] replaces: September 23, 2015 [Rev 2]