



NOBACERVIX[®]-sterile net wadding

REF 872406 (as of LOT 137166)

Product Description, Intended use, Application

The sterile, **cylindrical** tampon, **size 3.8 cm in diameter**, has a filling of absorbent, pressed cotton wool and a fine knit tubular wrapping. The extension of the knitted wrapping serves as a return strap. The tampon is suited for wound care in the case of gynaecological surgery. NOBACERVIX[®]-steril is for single use only and is individually sealed.

Composition

Cotton

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The tampon is tested according to DIN EN 14079 and USP/BP standards for cotton.

Sterilization complies with DIN EN ISO 11135.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: paper-film-packaging
Secondary packaging: folding box made of cellulose
Tertiary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.