



AUGENKLAPPE

REF 873010

Product Description, Intended use, Application

The rigid eye patch is used as protective bandage. The eye patch is fixed to the head by means of an elastic ribbon.

Composition

Cotton, cellulose, polyamide, elastomer

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE bag

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

Symbols used in labeling

Explanations at www.nobamed.com



Storage and transportation

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.