



## AUGENKOMPRESSE-steril

REF 873425

### Product Description, Intended use, Application

Individually packed sterile eye pad, oval shape, **53 x 70 mm**. The pad is die-cut from a cotton wool layer which is covered by gauze on both sides. Advantages of this light eye pad are its high air permeability as well as its softness. The eye pad supports eye dressing as absorbent, padding, protecting or caring part. AUGENKOMPRESSE-steril is for single use only.

### Composition

Cotton

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging:	paper-paper packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.