



## NOBALUMENAL®

REF 873550

### Product Description, Intended use, Application

The eye pad, size **54 x 70 mm** consists of a soft non-woven and is provided with an absorbent filling. It is closed all around so that no fluffs can emerge. Due to its **oval shape** the pad adapts to the eye socket. The eye pad supports eye dressing as absorbent, padding, protecting or caring part. NOBALUMENAL® is for single use only and must be sterilized with a validated procedure (radiation in accordance with DIN EN ISO 11137, or ethylene oxide in accordance with DIN EN ISO 11135).

### Composition

Polyester, viscose, rayon

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Sterilize before use on open wounds.

Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### Storage and Transport

Dry and dustfree, protect from sunlight.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.