



## NOBARHINAL<sup>®</sup>-basic

REF 874100

### Product Description, Intended use, Application

The nose bandage consists of a soft, cylindrical and absorbent padding (approx. 7.5 cm x 2 cm) that is covered by an elastic knit fabric. The elastic loops allow a fixation behind the ears. If used post-operatively, the nose bandage absorbs blood and wound secretions. It also serves as carrier material for aromatherapy.

### Composition

Polyamide, cotton

### Contra indication

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.