



NOBARHINAL[®]-basic large

REF 874110

Product Description, Intended use, Application

The nose bandage consists of a soft, cylindrical and absorbent padding (approx. 7.5 cm x 3 cm) that is covered by an elastic knit fabric. The elastic loops allow a fixation behind the ears. If used post-operatively, the nose bandage absorbs blood and wound secretions. It also serves as carrier material for aromatherapy.

Composition

Polyamide, cotton

Contra indication

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.