



NOBAMED Paul Danz AG

NOBARHINAL®

REF 874210

Product Description, Intended use, Application

The nose bandage (size L) consists of a soft, cylindrical and highly absorbent padding (approx. 8.5 cm x 3 cm) that is covered by an elastic knit fabric. The elastic loops allow a fixation behind the ears. If used post-operatively, the nose bandage absorbs blood and wound secretions. It also serves as carrier material for aromatherapy.

Composition

Polyamide, Cotton

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.