



## NOBA-ZUNGENLÄPPCHEN

REF 874800

### Product Description, Intended use, Application

NOBA® tongue pad, size **6 cm x 18 cm**, is made of 4-ply gauze. It is used to keep the tongue down in the case of laryngoscopy, thus preventing slipping during process. It is for single use.

### Composition

Cotton

### Notes

To increase patient and user safety, the tongue pad has been subjected to a validated germ reduction process and is available sterilized.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The dressing gauze complies with the normative requirements of DIN EN 14079 type 20.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging:	paper wrap
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.