



## NOBASTEP®

**REF 944600**

### **Product Description, Intended use, application**

The medical walking cast shoe (size XS, shoe size 31-34) is used with a thigh bandage or lower leg bandage (made of cast or of synthetic material). It enables walking with a support bandage and allows secure fixation thanks to the three Velcro fasteners NOBASTEP® may be used on the left or on the right. The material is water and dirt resistant. NOBASTEP® is washable at 30°C.

### **Composition**

PVC, TPR (thermoplastic styrene, butadiene plastic), polystyrole, naphthenic oil, calcium carbonate, zinc stearate, black dye), cardboard, EVA (ethylene vinyl acetate), polyester and nylon

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Carton made of cellulose

### **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

### **Storage and transport**

No special conditions

### **Single use device**

Intended for use on one patient. The re-processing of single use devices for reuse can decrease the product's performance significantly. Reusing a single use device can lead to microbiological danger.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.