



## WATTESTÄBCHEN

REF 974018

### Product Description, Intended use, Application

The cotton applicators (cotton tip: 10 - 11 mm) are used as applicator for liquids and ointments, for sanitation and for the collection of smear samples. For use on open wounds, the product must be sterilized with a validated sterilization method (Moist Heat in accordance with DIN EN ISO 17665, Radiation in accordance with DIN EN ISO 11137, or Ethylene Oxide in accordance with DIN EN ISO 11135) before. The cotton applicators have a rod made of polypropylene which is 15 cm long.

### Composition

Cotton, polypropylene

### Note

Sterilize before use on open wounds

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Legal Requirements

The cotton tipped applicator is a medical product according to the directive 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: PE Bag

Secondary packaging: Carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.