



## NOBATON®

REF 325050

### Product Description, Intended use, Application

The finely woven, tightly knitted tubular bandage, **size 5 cm x 25 m**, is made of unbleached cotton. A high crosswise elasticity ensures a good adaptability. The unbleached cotton shows water-repellent features and convinces by its excellent skin tolerance. NOBATON® serves as underlay under rigid or half-rigid cast bandages. Its tight structure protects the skin against material soaked with polyurethane resin and allows the application of a tightly fitting bandage. The rolled bandage is 25 m long. It is for single use and is individually packed.

### Composition

Cotton

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.