



NOBAMED Paul Danz AG

NOBATON®

REF 325120

Product Description, Intended use, Application

The finely woven, tightly knitted tubular bandage, **size 12 cm x 25 m**, is made of unbleached cotton. A high crosswise elasticity ensures a good adaptability. The unbleached cotton shows water-repellent features and convinces by its excellent skin tolerance. NOBATON® serves as underlay under rigid or half-rigid cast bandages. Its tight structure protects the skin against material soaked with polyurethane resin and allows the application of a tightly fitting bandage. The rolled bandage is 25 m long. It is for single use and is individually packed.

Composition

Cotton

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box
 made of
 cellulose

Secondary packaging: carton made of
 cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.