



## NOBATON<sup>®</sup>

### Product Description and Purpose

Finely woven, tightly knitted tubular bandage, which is made of unbleached cotton. A high crosswise elasticity ensures a good adaptability. The unbleached cotton shows water-repellent features and convinces by its excellent skin tolerance. NOBATON<sup>®</sup> serves as underlay under rigid or half-rigid cast bandages. Its tight structure protects the skin against the material soaked with polyurethane resin. It also allows the application of a tightly fitting bandage.

### Composition

Cotton

### Normative and Legal Requirements

NOBATON<sup>®</sup> is a medical product in accordance with the MDD, class I, rule 1, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### Storage

To be stored in a dry and dust-free environment

