



ABDOMINALVERBAND multisize

REF 371550

Product Description, Intended use, Application

The abdominal bandage, multisize 2: 95-125 cm, is permanently elastic and 25 cm wide. It is provided with a padding which is to lie on the abdomen. The hook and loop fastener adheres to the bandage on the full surface, so that there is a very large possibility of variation in the length. It serves to relieve the wound suture after operations in the abdominal cavity. The abdominal pad distributes the pressure on the abdominal wall evenly. The bandage is, fasteners closed, washable in the case of domestic use. Due to its material composition, the product may not be treated chemo-thermally and may therefore not be used on different patients. The abdominal bandage is not intended for permanent use as a dorsal bandage.

Composition

Cotton, polyamide, viscose, polyester, elastane

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: box made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

Storage and Transport

Dry and dustfree.

Single use device

Intended for use on one patient. The re-processing of single use devices for reuse can decrease the product's performance significantly. Reusing a single use device can lead to microbiological danger.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.