



NOBALUX®

REF 416365

Product Description, Intended use, Application

This white bandage with a **width of 6 cm** and a **stretched length of 4 m** is an elastic dressing retention bandage. Due to interlaced yarns, the fixation bandage is lengthwise elastic as well as crosswise elastic. NOBALUX® is especially suited for cases where the fixation of the wound pad is extremely difficult. NOBALUX® is for single use only.

Composition

Viscose, polyamide

Contra-Indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.