



## NOBASTRETCH®-kräftig

REF 481708

### Product Description and Purpose

The brown long-stretch compression bandage (size: 8 cm x 7 m) has a stretchability of approx. 170 % and is permanently elastic. Owing to its weave, it has a strong structure. Long-stretch bandages are mainly used to create bandages that will support and relieve muscles or joints in cases of disorders, for example in the case of sports injuries. They are also used in cases of phlebologic indications like superficial varicose veins and for the follow-up treatment of venous conditions. Long-stretch bandages may also be used to complement compression bandages during the day to flush out oedema. If so used, and if used as phlebologic compression bandage in general, the relevant medical guidelines must be observed. Due to the high resting pressure NOBASTRETCH®-kräftig is not left applied overnight or in the case of longer periods of rest during the day.

### Contra-Indications

Not be used if there is a known allergy to one of the below-mentioned components. Do not use in the case of deep venous disorders. In the case of a medical compression therapy the contra-indications and risks mentioned in the relevant medical guidelines must be observed.

### Composition

Cotton, polyamide, polyurethane

### Normative and Legal Requirements

NOBASTRETCH®-kräftig is a medical device according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### Packaging

Primary packaging: folding box  
made of cellulose

Secondary packaging: carton  
made of cellulose

### Storage

To be stored in a dry and dust-free environment

### The product bears following labels and symbols:

