



NOBACREPP®

REF 490403

Product Description, Intended use, Application

The textile elastic fast edge bandage (width: 2.5 cm) consists of unbleached cotton and is used for fixation and supporting. Due to the way it is manufactured, it has a terrycloth like fabric structure. This makes it particularly suitable for the surgical field (e.g. hand surgery), as it has a padding effect. The bandage is 4 m long and is individually packed.

Composition

Cotton

Contraindications

The product should not be used in the case of a known allergy against the material.

Note

The product has been subject to a validated microbial reduction process and is therefore pre-sterilized

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

Packaging

Primary:	folding box made of cellulose
Secondary:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.