



## NOBACREPP®

REF 490408

### Product Description, Intended use, Application

The textile elastic fast edge bandage (width: 8 cm) consists of unbleached cotton and is used for fixation and supporting. Due to the way it is manufactured, it has a terrycloth like fabric structure. This makes it particularly suitable for the surgical field (e.g. hand surgery), as it has a padding effect. The bandage is 4 m long and is individually packed.

### Composition

Cotton

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Note

The product has been subject to a validated microbial reduction process and is therefore pre-sterilized

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

### Packaging

Primary:	folding box made of cellulose
Secondary:	carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.