



NOBAPAD®

REF 600804

Product Description and Purpose

Especially soft and stretchable padding bandage, **size 4 cm x 3 m**, which is made of binder-free synthetic wadding. NOBAPAD® can be used as preparatory padding for treatments with plaster, polyester or fibreglass support bandages. It is likewise suited for the padding of Cramer splints, Volkmann splints or other unpadded splints, or else as padding underneath compression bandages and as padding material for tubular bandages. The rolled bandage is for single use and is **individually packed in a paper banderole**.

Contra-Indications

No contra-indications known.

The product should not be used if there is a known hypersensitivity to the below-mentioned ingredient.

Composition

Polyester

Normative and Legal Requirements

NOBAPAD® is a class I, rule 1 medical device according to the MDD, the Directive 93/42/EEC, and the Regulation MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: paper

Secondary packaging: folding box
made of cellulose

Tertiary packaging: carton
made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

