



## NOBAPAD®

REF 600820

### Product Description and Purpose

Especially soft and stretchable padding bandage, **size 20 cm x 3 m**, which is made of binder-free synthetic wadding. NOBAPAD® can be used as preparatory padding for treatments with plaster, polyester or fibreglass support bandages. It is likewise suited for the padding of Cramer splints, Volkmann splints or other unpadded splints, or else as padding underneath compression bandages and as padding material for tubular bandages. The rolled bandage is for single use and is **individually packed in a paper banderole**.

### Contra-Indications

No contra-indications known.

The product should not be used if there is a known hypersensitivity to the below-mentioned ingredient.

### Composition

Polyester

### Normative and Legal Requirements

NOBAPAD® is a class I, rule 1 medical device according to the MDD, the Directive 93/42/EEC, and the Regulation MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### Packaging

Primary packaging: paper

Secondary packaging: folding box  
made of cellulose

Tertiary packaging: carton  
made of cellulose

### Storage

To be stored in a dry and dust-free environment

**The product bears following labels and symbols:**

