



## UNTERSUCHUNGSKITTEL

REF 628105

### Product Description, Intended Use, Application

Disposable examination gown (size M) for use in medical and nursing activities, mainly in ambulances to protect against exposure to infectious agents. It is also intended to protect from a potential exposure to small quantities of spray or accidental low volume splashes of liquid chemicals.

The gown is closed in the back and is equipped with cuffs.

### Composition

Polypropylene, polyethylene, polyester

### Note

The product has been subjected to a validated microbial reduction process and is therefore pre-sterilized

### Normative and Regulatory, Requirements, Common Standards

The gown complies with the European standards for Personal Protective Equipment in accordance with the provisions of the EU Regulation 2016/425 category III-module C2.

The gown protects against infections in accordance with EN ISO 14126:2003 + AC:2004 and protects against liquid chemical, light spray (type 6B) according to EN 13034:2005 + A1:2009

EN ISO 13688: 2013+ A1: 2021

EN 14325: 2004

Further information can be found in the instruction for use.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging:

OPP pouch

Secondary packaging:

carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



CE 0624 PSA (CAT III)

### Storage and Transport

Dry and dustfree, away from sources of light and heat

### Disposal

If not contaminated the product can be treated as a common textile waste. If contaminated, it should be treated as harmful garbage and discarded according to country laws.