

NOBAMED Paul Danz AG

UNTERSUCHUNGSKITTEL

REF 628120

Product Description and Purpose

Disposable examination gown (size XXL) for use in medical and nursing activities, mainly in ambulances to protect against infectious exposure to agents. The product is liquid tight. The gown is closed in the back and is equipped with cuffs.

Composition

Polypropylene, Polyethylene

Normative and Legal Requirements

Examination gowns are protective equipment according to the PPE Regulation (EU) 2016/425 category III.

The gown protects against infections in accordance with

EN ISO 14126: 2003 + AC: 2004,

ISO 16603/16604,

ISO 22610,

ISO 22611 and

ISO 22612.

Further information can be found in the product information.

The product does not contain dangerous toxic substances according to REACH.

All packing levels are labelled acc. to DIN EN ISO 15223-1 and EN 1041.

Packaging

Primary packaging: OPP pouch

Secondary packaging: carton made

Of cellulose

Storage

To be stored in a dry and dust-free environment.

The product bears the following symbols and marking







C € 0624 PPE (CAT III)

Date of information: May 12, 2020 [REV 0]