



NOBAMED Paul Danz AG

NOBATRAST® VAKUUM-OP-SET WUNDHAKENBEZUG

REF 644624

Product Description and Purpose

Sterile disposable single-layer white covers sewn of medically bleached knitted cotton fabric, 4.5 cm x 24.0 cm, packed per 6 pcs, with worked-in X-ray thread and with bands. The bands are used to fix the covers at the wound hook. The wound hook covers are put directly on the sterile wound hooks. They reduce light reflections caused by the metal due to the operation lamp. Sensitive tissue is protected against the instruments' sharp edges. The covers also protect the instruments until they are used on the patient. The wound hook covers are vacuum-packed in two transparent film-film-packages. They are firmly positioned in the package. Owing to a counting card on the inner packaging and a second counting card on the inner side of the outer packaging, the counting may be checked quickly and safely. The vacuum of the inner and outer packaging guarantees the sterility of the set until it is used.

Contraindication:

Known allergy against cotton

Composition

Cotton,
X-Ray thread contains barium sulphate,
PVC and DINP

Normative and Legal Requirements

NOBATRAST® VAKUUM-OP-SET
WUNDHAKENBEZUG is a medical product according to the MDD 93/42/EEC and the MDR (EU) 2017/745 and is classified as class IIa, rule 7.

Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: film-film packaging

Secondary packaging: film-film packaging

Tertiary packaging: folding box made of cellulose

Quaternary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

The product bears the following labels and symbols:



STERILE R



CE
0123