



## NOBATRAST<sup>®</sup>-Vakuum-OP-Set Präpariertupfer

REF 645610

### Product Description and Purpose

The double packed preparation balls made of 24-thread dressing gauze, 10 cm x 10 cm, (size L), packed per 10 pcs, have a woven-in X-ray thread. They are used to prepare delicate tissue in the surgical field (for example blood vessels, fine tendons). The balls are separately placed in a transparent tray, thus enabling a visual content check before the operation. Owing to a yellow counting card on the inner packaging and a second yellow counting card on the inner side of the outer packaging, the counting may be checked quickly and safely. The vacuum of the inner and outer packaging guarantee the sterility of the set until it is used.

### Contraindication

Allergy against cotton, barium sulfate, PVC or DINP

### Composition

Cotton, X-Ray thread contains barium sulfate, PVC and DINP

### Normative and Legal Requirements

NOBATRAST<sup>®</sup> Vakuum-OP-Set Präpariertupfer (Vakuum OP Set preparation balls) is a medical product according to the MDD 93/42/EEC and the MDR (EU) 2017/745 and is classified as class IIa, rule 7.

The gauze complies with the requirements of DIN EN 14079 type 24.

Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### Packaging

Primary packaging:	film-film packaging
Secondary packaging:	film-film packaging
Tertiary packaging:	folding box made of cellulose
Quaternary packaging:	carton made of cellulose

### Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

The product bears the following symbols and marking:

