



NOBAMED Paul Danz AG

NOBATRAST® VAKUUM-OP-SET NOBAOPHTHALM®

REF 647110

Product Description and Purpose

NOBATRAST® VAKUUM-OP-SET NOBAOPHTHALM® contains ophthalmological arrow sponges with bevel cut that are used for preparation and as swabs in ophthalmology. They have a high suction speed and excellent absorbency. The sponge swells up when liquids are adsorbed and retains them.

The set is supplied in double sterile packaging.

Composition

Viscose nonwoven, cellulose

Normative and Legal Requirements

NOBATRAST® VAKUUM-OP-SET NOBAOPHTHALM® is a medical product in accordance with the MDD, class IIa, rule 6, and complies with the requirements of the directive 93/42/EEC.

Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: film-film packaging

Secondary packaging : film-film packaging

Tertiary packaging: folding box made of cellulose

Quaternary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

The product bears the following labels and symbols

