



NOBAMED Paul Danz AG

WASCHHANDSCHUHE

laminated

REF 672300

Product Description, Intended use, Application

The soft outside of the washing gloves consists of polyester, whereas the inner side is plasticized to prevent the contact with the washing water. The washing gloves are gentle to the patient's skin, and the nursing person is protected against water and against the patient's skin germs. The gloves are disposable and suited for washing MRSA patients if examination gloves are also worn. The WASCHHANDSCHUH is used to clean the patient's body to prevent the spread of germs.

Composition

Polyester, polyethylene

Contra-Indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.