



## **BARETTHAUBE basic green**

**REF 674110**

### **Product Description and Purpose**

The caps are made of a light non-woven which is air-permeable and breathable. In this way a heat accumulation is prevented and the wearing comfort is increased. The bouffant caps cover the hair completely and reliably and retain skin and hair particles.

### **Composition**

Polypropylene, Elastomer thread (latex-free)

### **Normative and Legal Requirements**

Baretthaube (bouffant caps) are medical products according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and are classified as class I, rule 1

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### **Packaging**

Primary packaging: folding box made of Cellulose

Secondary packaging: carton made of cellulose

### **Storage**

To be stored in a dry and dust-free environment

**The product bears following labels and symbols:**

