



NOBAMED Paul Danz AG

UNIVERSALKLEMME

REF 682500

Product Description, Intended use, Application

Anatomically shaped disposable swab clamp/tube clamp, 12 cm, made of blue plastic, with locking mechanism. The clamp is for universal use in patient care in hospital wards, out-patient facilities, or medical practices. The serrated part serves as a swab holder. The rear part is suitable for catheter clamping.

Composition

Polypropylene

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary	PE bag
Secondary:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and ISO 20417.

Storage and transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.