

NOBAMED Paul Danz AG

NOBAPEG®

Product Description and Purpose

NOBAPEG® is a disposable medical device combination with two non-woven slit swabs (NOBATOP®-S) 7.5 x 7.5 cm, three gauze swabs 7.5 x 7.5 cm, 8-ply, as well as an adhesive non-woven 12 x 15 cm. The product is perfectly suited for an appropriate change of dressing at the point of exit of the tube and place of entry of the catheter, in particular for percutaneous endoscopic gastrostomy (PEG) tubes for enteral feeding.

Composition

Cotton, viscose, polyester, polyethylene, polyacrylate adhesive

Normative and Legal Requirements

NOBAPEG® is a medical product in accordance with the MDD, class Is, rule 4, and complies with the requirements of the directive 93/42/EEC.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: paper-film

packaging

Secondary packaging: folding box made

of cellulose

Tertiary packaging: carton made of

cellulose

Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

The product is bearing the C € marking, in compliance with the directive's provisions











Date of information: September 25, 2015 [REV 3]

replaced: March 4, 2015