



VERBANDWECHSEL-SET

REF 696000 (as of LOT 137704)

Product Description, Intended use, Application

The sterile VERBANDWECHSEL-SET consists of 5 binder-free, 4-ply non-woven swabs NOBATOP® 8, size 10 x 10 cm and size 5 x 5 cm. The swabs are packed in a white plastic bowl, which enables a sterile change of dressing. The non-woven swabs are used in general wound care, for wound cleansing, wound covering and for the absorption of blood and wound fluids.

Composition

NOBATOP® 8: viscose, polyester
Plastic bowl: PET

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

Sterilization complies with DIN EN ISO 11135
The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging:	paper-film packaging
Secondary packaging:	box made of cellulose
Tertiary packaging	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417.

Storage and Transport

Dry and dustfree

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.