



MULLTUPFER-steril

REF 735725

Product Description and Purpose

The sterile gauze balls have a diameter of 3.5 cm. The gauze (DIN EN 14079 type 20) is twisted into a round shape without using any rubber ring. They are used for the treatment of superficial wounds and for skin disinfection. One sterile pouch contains 20 gauze balls.

Contra-Indications

No contra-indications known. The product should not be used in case of a known allergy against the component mentioned below.

Composition

Cotton

Normative and Legal Requirements

Mulltupfer-steril (sterile gauze balls) are medical products according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and are classified as class Is, rule 4.

The gauze complies with the requirements of DIN EN 14079, type 20.

Sterilization complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	film-film packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	Carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

