



NOBAMED Paul Danz AG

NOBARAPID®

REF 773107

Product Description and Purpose

Wound swab, size **7.5 cm x 7.5 cm**, individual sterile packaging. 50 sterile pouches are packed in a folding box. NOBARAPID® wound swabs consist of an absorbent core made of cotton-polyester which has been laminated with a perforated polyethylene film on both sides. The films prevent the swabs from adhering to the wound. NOBARAPID® wound swab is suitable for low to moderately exuding wounds, post-operative wounds and painful wounds like skin abrasions and cuts, for example.

They are for single use only.

Contra-Indications

No contra-indications known.

The product should not be used if there is a known hypersensitivity to one of its ingredients.

Composition

Cotton, polyester and polyethylene

Normative and Legal Requirements

NOBARAPID® is a class Is, rule 4 medical device according to the MDD, the directive 93/42/EEC and the MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN ISO 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: paper packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

