



NOBAMED Paul Danz AG

NOBACOLLOID®

REF 790105

Product Description and Purpose

Individually packed, sterile hydroactive wound dressing for moist wound treatment, 5 cm wide and 10 cm long. The hydrocolloid matrix provides a very high absorbency of the wound exudate, causing the formation of a cohesive gel, while the breathable polyurethane film allows an excellent moisture vapor transmission. In this way an optimal moist wound environment is ensured which supports the wound healing process, and the risk of damaging the intact periwound area is minimized. NOBACOLLOID® can be left in place for up to seven days. The wound dressing may be used for the management of light to moderately exuding wounds, such as pressure sores and leg ulcers, as well as partial and full thickness wounds, donor sites, abrasions, trauma wounds, post-operative surgical wounds and superficial burns (1st and 2nd degree).

Composition

Hydrocolloid matrix, polyurethane film

Normative and Legal Requirements

NOBACOLLOID® is a medical product according to the MDD 93/42/EEC and the MDR (EU) 2017/745 and is classified as class IIb, rule 4. Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: paper-foil packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose


Storage

To be stored in a dry and dust-free environment

The product bears the following symbols and marking:



STERILE R



CE

0123