



NOBAMED Paul Danz AG

NOBAMED® TUPFER

REF 880015

Product Description and Purpose

The unsterile NOBAMED-balls, 3.0 cm diameter, are made of soft, absorbent, binder-free produced non-woven fabric and contain a silicone ring. They are used for wound care and in combination with a disinfectant for the disinfection of the skin. Prior to these activities, the product has to be sterilized with a validated procedure (moist heat at 121°C / 134°C, 2 to 3 bar), according to DIN EN ISO 17665, radiation according to DIN EN ISO 11137 or ethylene oxide according to DIN EN ISO 11135)

Composition

Viscose, polyester, silicone

Normative and Legal Requirements

NOBAMED® TUPFER is a medical product according to the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is, according to rule 4.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:



Date of information: April 01,2020 [REV4]
replaced: September 25, 2015 [REV3]