



**NOBAMED Paul Danz AG**

## **NOBAMED® TUPFER**

**REF 880020**

### **Product Description and Purpose**

The unsterile NOBAMED-balls, 3.5 cm diameter, are made of soft, absorbent, binder-free produced non-woven fabric and contain a silicone ring. They are used for wound care and in combination with a disinfectant for the disinfection of the skin. Prior to these activities, the product has to be sterilized with a validated procedure (moist heat at 121°C / 134°C, 2 to 3 bar), according to DIN EN ISO 17665, radiation according to DIN EN ISO 11137 or ethylene oxide according to DIN EN ISO 11135)

### **Composition**

Viscose, polyester, silicone

### **Normative and Legal Requirements**

NOBAMED® TUPFER is a medical product according to the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is, according to rule 4.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### **Packaging**

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

### **Storage**

To be stored in a dry and dust-free environment

**The product bears following labels and symbols:**



Date of information: April 01,2020 [REV4]  
replaced: September 25, 2015 [REV3]