



ARMTRAGETUCH

REF 900150

Product Description and Purpose

The triangular bandage is made of black viscose fabric and is used for the short term immobilization of the forearm and for the fixation of wound pads in a first aid setting.

Composition

Viscose

Normative and Legal Requirements

The triangular bandage is a medical product according to the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, rule 1

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PE pouch

Secondary packaging: Carton
made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

