



## Armtragetuch

**REF 900150**

### Product Description, Intended use, Application

The triangular bandage, 96 cm x 96 cm x 136 cm, is made of white viscose fabric and is used for the short term immobilization of the forearm and for the fixation of wound pads in a first aid setting.

### Composition

Viscose

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDR (EU) 2017/745*.

The bandage complies with DIN standard 13168.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary PE-pouch

Secondary: carton  
made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.