



ARMTRAGETUCH

REF 900200

Product Description and Purpose

The triangular bandage, 96 cm x 96 cm x 136 cm, is made of white nonwoven fabric and is used for the short term immobilization of the forearm and for the fixation of wound pads in a first aid setting.

Composition

Polypropylene

Normative and Legal Requirements

The triangular bandage is a medical product according to the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, rule 1

The bandage complies with DIN standard 13168.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PE pouch

Secondary packaging: Carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

