



## Armtragegurt

REF 910035

### Product Description, Intended use, Application

The black arm sling support (width: 35 mm) is provided with a double sling and a clasp which is adjustable. By means of the product the forearm can be immobilized for a short time.

### Composition

Cotton, stainless steel

### Contraindication

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

|            |                          |
|------------|--------------------------|
| Primary:   | PE bag                   |
| Secondary: | carton made of cellulose |

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Single use device

Intended for use on one patient. The re-processing of single use devices for reuse can decrease the product's performance significantly. Reusing a single use device can lead to microbiological danger.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.