



NOBAMED Paul Danz AG

MUNDSPATEL-steril

REF 922600

Product Description and Purpose

The tongue depressor (size 18 mm x 150 mm) is made of clean processed wood which is smoothly polished. It serves to inspect the pharyngeal space and to apply ointment. Due to its sterility, the tongue depressor, which is individually packed, can also be used under aseptic conditions.

Contraindications

Known allergy against the component listed below

Composition

Wood

Normative and Legal Requirements

The tongue depressor is a medical product according to the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is, according to rule 5.

Sterilization of the product is in accordance with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH. It has DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	paper/ foil packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears the following symbols and marking:

