



WATTESTÄBCHEN-steril

REF 974411

Product Description and Purpose

The sterile cotton applicators (cotton tip: 10 - 11 mm) are used to apply liquids (for example wound disinfectants) and ointments, to clean the wound by means of a suitable solution and for collection of smear samples. The cotton applicator has a rod made of polypropylene which is 15 cm long. One sterile pouch contains 2 cotton applicators.

Contraindications

Not to be used in case of a known allergy against one of the components listed below

Composition

Cotton, polypropylene

Normative and Legal Requirements

WATTESTÄBCHEN-steril is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is, rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	paper- film packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment.

The product bears the following labels and symbols:

